



RISK CALCULATION SHOULD TAKE LESS THAN 9 MONTHS

– P5Link assists you with saving time –

Typolog
software

P5LINK

P5LINK

P5Link allows fully automated risk calculation with Prisca 5®, the prenatal screening software, distributed by Siemens Healthineers.

P5Link is **compatible with each Laboratory Information System** (LIS) which is capable of generating ASTM files. Importing all data into Prisca, using Prisca to calculate risks and export all results back to the LIS **takes a few seconds only.**

Thus, P5Link

- » saves valuable time
- » helps avoiding typing errors
- » can process twin and egg donor pregnancies

Three modes are available:

- » Full automation with validation of the results by P5Link
- » Full automation with validation of the results by the customer
- » Semiautomatic mode

P5Link is the latest software developed by Typolog Software. Based on the current scientific status and **experience of more than 20 years in prenatal diagnosis.**



P5LINKPREPARE

P5LinkPrepare is a configuration tool which guides you through the process of setting up P5Link in such a way, that it is able to communicate with your LIS and with Prisca.

With P5LinkPrepare these steps **become very easy and straightforward** and can be done with a minimum of time.

Minimum requirements for P5Link:

- » Windows 7 or higher
- » Screen resolution 1.440 * 900 pixels or higher
- » 500 MB free hard disk space
- » Prisca version 5.X



TYOLOG SOFTWARE

Typolog Software cooperates with **an international team, comprising mathematicians, biologists and IT experts.** Our key expertise is in designing, managing and developing software, specialised in the medical and statistical area.

Having worked in the area of medical informatics for more than 20 years we are **internationally renowned for our expertise and excellent service.**

Ever since the launch of PRISCA – the world-wide first Windows program for the Triple-Test – we have been developing tools, which **suit our customers' needs - with „creativity and integrity“.**

Typolog Software Ltd. & Co KG
Esinger Straße 69B
25436 Tornesch
Germany
Fax: +49-4122-960 436
Email: medical@typolog.de